

**APPLICATION NOTE** FOR THERMO SCIENTIFIC QE ORBITRAP™ MS

# Localization of Multiple Phosphorylation Sites in the case of $\alpha$ S-Casein

Yury V. Vasil'ev and Jared B. Shaw

## Introduction

Casein phosphorylation is a post-translational modification (PTM) catalyzed by kinase enzymes when phosphate groups are attached to specific amino acids<sup>1</sup>. The labile nature of phosphorylation modifications makes it very difficult to study them with traditional tandem mass spectrometry (MS/MS) techniques using collision-induced dissociation (CID). Electron capture dissociation (ECD) is known as a soft MS/MS technique with ability to preserve labile PTM during fragmentation. The present work shows using ECD method for localization of phosphorylation in the case of  $\alpha$ S-casein.


## Experimental

### Sample preparation

Phosphorylated  $\alpha$ S-casein from bovine milk was bought from Sigma-Aldrich and used without further purification. The  $\alpha$ S-casein sample of 1  $\mu$ M concentration was dissolved in a 1:1 (v:v) water-methanol mixture with 0.1% concentration of formic acid.

### Instrumentation


ECD spectra of the  $\alpha$ S-casein were obtained with an ExD cell (e-MSion, Inc.)<sup>2</sup> installed on a Thermo Scientific Q Exactive Orbitrap mass spectrometer (Figure 1). The sample was introduced through direct infusion at a flow rate of 5  $\mu$ L per minute. The autotune feature in the ExDControl software (e-MSion, Inc.) was used to optimize the ExD cell voltage profiles for transmission and for ECD fragmentation, both with the filament heating current set to 2.2 A. To tune for transmission, Thermo-Fisher calibrant solution was used, whereas tuning for ECD the peptide standard Substance P and bovine ubiquitin were used. To preserve phosphorylation in the  $\alpha$ -casein sample, ECD spectra with low HCD supplemental activation (1 and 5 V) were recorded. To analyze the EChcD spectra, Viewer software (e-MSion, Inc.) was used at 10 ppm tolerance for masses.




**Figure 1.** Schematic of the ExD cell in the Thermo Scientific Q Exactive Orbitrap mass spectrometer.

## Results and Discussion


Figure 2 shows the ECD spectrum of  $\alpha$ -casein obtained with HCD supplemental activation 5 V. All eighth phosphate groups attached to some serine residues were identified in the spectrum (see some of them in Figure 3) with total sequence coverage 95% as was determined by Viewer (Figure 4).



**Figure 2.** EChcD spectrum of  $\alpha$ -casein for 20+ charge state precursor. Blue filled circles show ECD-type *c*- and *z*-fragments, green filled circles are CID-type *b*- and *γ*-fragments, purple filled circles are *a*-ions, open circles are secondary *d*- and *w*-ions.



**Figure 3.** Some ECD-fragments with phosphoserines:  $c_{46}^{6+}$  (A),  $c_{48}^{6+}$  (B),  $c_{67}^{5+}$  (C),  $c_{75}^{6+}$  (D).



**Figure 4.** Sequence coverage determined by Viewer from EChcD spectrum of  $\alpha$ S-casein in Figure 2. Serines with phosphate groups are shown in red.

## References

1. Bijl, E., van Valenberg, H.J.F., Huppertz, T., van Hooijdonk, A.C.M., Bovenhuis, H. (2014) Phosphorylation of  $\alpha$ S1-casein is regulated by different genes, *J. Dairy Sci.* 97 (11), 7240-6.
2. Shaw, J. B.; Malhan, N.; Vasil'ev, Y. V.; Lopez, N. I.; Makarov, A.; Beckman, J. S.; Voinov, V. G. Sequencing Grade Tandem Mass Spectrometry for Top-Down Proteomics Using Hybrid Electron Capture Dissociation Methods in a Benchtop Orbitrap Mass Spectrometer. *Anal. Chem.* **2018**, 90 (18), 10819-10827. <https://doi.org/10.1021/acs.analchem.8b01901>.

### SALES

---

e-MSion, Inc.

2121 NE Jack London Street

Suite 140

Corvallis, Oregon 97330

[Sales@e-MSion.com](mailto:Sales@e-MSion.com)



<https://e-msion.com/products/thermo-fisher-scientific/>

